

Suction Irrigation Device

Stapleline[®]
Quality that connects



CE 0481

5 mm laparoscopic Suction Irrigation Device

- the instrument fits the surgeon's palm properly
- the same operating comfort for right and left handers
- suction and irrigation combination easy to handle with an action radius of 360°
- enormous throughput by larged-bore valve system in the handle
- the tube is 3 m in length

Instrument No.	Size	Sales Unit
CX-01	5mm x 230mm	10
CX-02	5mm x 330mm	10
CX-03	5mm x 430mm	10

Stapleline Medizintechnik GmbH

+49 (0) 234 - 936 487 10

+49 (0) 234 - 936 487 12

Bessemer Straße 30
44793 Bochum, Germany 

www.stapleline.com

E-Mail: info@stapleline.com