

Skin Stapler & Remover

CE 0197



Single-use Skin Stapler for wound closure and surgical incisions. Using the Skin Stapler Remover the skin staples can be easily removed.

Material

The staples consist of stainless steel according to AISI 316L.
The handhold is made of ABS according to CAS 9003-56-9. The Skin Stapler Remover consists of ABS and stainless steel 1Cr18Ni9.

Performance features

Sterile, excellent biocompatibility, very good metal stability, quick realization of the surgical suture, clear view to the wound edges, suitable for all cross sections.

Instrument No.	Size of Staples	Sales Unit
Skin Stapler: 35W	6.4mm x 4.0mm	100
Remover: SRK-T-2		200

Stapleline Medizintechnik GmbH

+49 (0) 234 - 936 487 10

+49 (0) 234 - 936 487 12

Bessemer Straße 30
44793 Bochum, Germany

www.stapleline.com

E-Mail: info@stapleline.com